

## The 3-Splint-System

Three material thicknesses, worn during 4 weeks for a unique and beautiful smile: With each treatment step, you'll be wearing 3 splints of different thicknesses.

Thanks to this steady increase in pressure, your teeth are aligned efficiently, yet smoothly and with extraordinary comfort.

For easy distinction, the corresponding degree of hardness is marked on the aligners in the posterior segment.

We look forward to making you smile!



## Introducing the CA<sup>®</sup> CLEAR ALIGNER



STEP 1	Week	Material Thickness
STEP 1	Week 1	soft
	Week 2	medium
	Week 3	hard
	Week 4	hard
STEP 2	Week 5	soft
	Week 6	medium 0.625 mm
	Week 7	hard
	Week 8	hard
STEP 3	Week 9	soft
	Week 10	medium 0.75 mm
	Week 11	hard
	Week 12	hard



Suite 2B, 1-3 Gurrigal St  
Mosman, NSW, 2088

Level 5, 13 Spring St  
Chatswood, NSW, 2067

Ph: 1300 354 665  
Email: [reception@orthodontistsydney.com](mailto:reception@orthodontistsydney.com)  
Website: [www.orthodontistsydney.com](http://www.orthodontistsydney.com)

### The latest in tooth aligning technology

- **Transparent**  
Virtually invisible for easy wear
- **Customised**  
Individually planned and constructed by orthodontists and orthodontic laboratories
- **Comfortable**  
Smooth finish for optimum comfort
- **Removable**  
Remove and re-insert at your leisure
- **Effective**  
Precise movements calculated by orthodontists for a beautiful smile

Image courtesy of SCHEUDENTAL GmbH

## What is CA<sup>®</sup> CLEAR ALIGNER?

A series of clear, custom made, removable splints which uses gentle forces to align teeth to the desired position. The CA<sup>®</sup> material is fabricated to produce lower forces to ensure a gentle yet fast tooth alignment. The CA<sup>®</sup> CLEAR ALIGNER can be worn with confidence, patients can be assured that when wearing the splint they can speak with ease, and best of all feel confident when smiling and engaging in everyday activities.



CA<sup>®</sup> CLEAR ALIGNER splints are virtually invisible and so easy to insert and remove.



## How does CA<sup>®</sup> CLEAR ALIGNER work?

Each CA<sup>®</sup> splint is custom made to ensure comfort with specific and precise tooth movements calculated by Dr. Jones using the latest technology. All CA<sup>®</sup> CLEAR ALIGNERS are fabricated in Dr. Jones' lab on site which allows him to make necessary adjustments and to have a prompt turn-around time for his patients.

The CA<sup>®</sup> splints are to be worn for a minimum of 22 hours per day to ensure optimum results. Due to the gentle and comfortable nature of the tooth alignment this can be easily achieved by patients.

Patients undergoing treatment with CA<sup>®</sup> CLEAR ALIGNER are monitored closely by Dr. Jones, with an impression being taken after each set of aligners to ensure optimum tooth movement. With state of the art technology Dr. Jones is able to make a video of the calculated tooth movement and the determined amount of CA<sup>®</sup> splints needed to reach the desired tooth correction, he can then share this with his patients to ensure communication is open and treatment goals are fully understood by the patient.



## Dr. Jones and the CA<sup>®</sup> system

Dr. Jones understands that conventional methods of tooth alignment are not feasible for some patients, which is why he is excited to be initiating the CA<sup>®</sup> CLEAR ALIGNER system.

Dr. Jones has a long history with aligners and his expertise in this field has allowed him to travel the globe, constantly learning and collaborating with like-minded professionals to structure innovative methods of tooth alignment to integrate into his high quality orthodontic treatments. With recent travel to Germany to extend his knowledge on the CA<sup>®</sup> system Dr. Jones found that with this latest technology he has greater control over tooth alignment than he ever did with other removable tooth straightening systems, with more flexibility and greater possibilities for his patients.

## Treatment with CA<sup>®</sup> CLEAR ALIGNER

The CA<sup>®</sup> system is mostly used for minor tooth movement or relapse, however, it can also be used in conjunction with braces to minimise time spent in braces. The more tooth movement required, the more CA<sup>®</sup> CLEAR ALIGNER splints will be required.

Each patient's treatment is individualised so it is best to make an appointment with Dr. Jones to discuss your treatment goals.

Find out more at:  
[www.ca-clear-aligner.com](http://www.ca-clear-aligner.com)